



PCOLCE Antibody, HRP conjugated

Product Code	CSB-PA614881LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15113
Immunogen	Recombinant Human Procollagen C-endopeptidase enhancer 1 protein (169-449AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Procollagen C-endopeptidase enhancer 1 (Procollagen COOH-terminal proteinase enhancer 1) (PCPE-1) (Procollagen C-proteinase enhancer 1) (Type 1 procollagen C-proteinase enhancer protein) (Type I procollagen COOH-terminal proteinase enhancer), PCOLCE, PCPE1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PCOLCE
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.