



GRM3 Antibody, Biotin conjugated

Product Code	CSB-PA622790LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14832
Immunogen	Recombinant Human Metabotropic glutamate receptor 3 protein (46-354AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Metabotropic glutamate receptor 3 (mGluR3), GRM3, GPRC1C MGLUR3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GRM3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.