



# CARD16 Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA682727LB01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q5EG05
<b>Immunogen</b>	Recombinant Human Caspase recruitment domain-containing protein 16 protein (1-197AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Caspase recruitment domain-containing protein 16 (Caspase recruitment domain-only protein 1) (CARD-only protein 1) (Caspase-1 inhibitor COP) (Pseudo interleukin-1 beta converting enzyme) (Pseudo-ICE) (Pseudo-IL1B-converting enzyme), CARD16, COP COP1
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cell Biology
<b>Target Names</b>	CARD16