



omp Antibody, HRP conjugated

Product Code	CSB-PA706583LB01RMU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q52764
Immunogen	Recombinant Rickettsia japonica 17 kDa surface antigen protein (20-159AA)
Raised In	Rabbit
Species Reactivity	Rickettsia japonica
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	17 kDa surface antigen, omp
Immunogen Species	Rickettsia japonica (strain ATCC VR-1363 / YH)
Research Area	Microbiology
Target Names	omp
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.