



KCNV1 Antibody, HRP conjugated

Product Code	CSB-PA747493LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6PIU1
Immunogen	Recombinant Human Potassium voltage-gated channel subfamily V member 1 protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Potassium voltage-gated channel subfamily V member 1 (Neuronal potassium channel alpha subunit HNKA) (Voltage-gated potassium channel subunit Kv8.1), KCNV1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KCNV1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.