



DNER Antibody, HRP conjugated

Product Code	CSB-PA818784NB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NFT8
Immunogen	Peptide sequence from Human Delta and Notch-like epidermal growth factor-related receptor protein (150-168AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Delta and Notch-like epidermal growth factor-related receptor, DNER, BET
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	DNER
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.