



# APH1B Antibody, Biotin conjugated

Product Code	CSB-PA819900ND01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WW43
Immunogen	Peptide sequence from Human Gamma-secretase subunit APH-1B protein (93-110AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gamma-secretase subunit APH-1B (APH-1b) (Aph-1beta) (Presenilin-stabilization factor-like), APH1B, PSFL
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	APH1B