



RORB Antibody, Biotin conjugated

Product Code	CSB-PA821694LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92753
Immunogen	Recombinant Human Nuclear receptor ROR-beta protein (1-269AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear receptor ROR-beta (Nuclear receptor RZR-beta) (Nuclear receptor subfamily 1 group F member 2) (Retinoid-related orphan receptor-beta), RORB, NR1F2 RZRB
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RORB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.