



NCLN Antibody, Biotin conjugated

Product Code	CSB-PA822169LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q969V3
Immunogen	Recombinant Human Nicalin protein (214-522AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nicalin (Nicastrin-like protein), NCLN
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NCLN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.