



HTRA1 Antibody, Biotin conjugated

Product Code	CSB-PA835695LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92743
Immunogen	Recombinant Human Serine protease HTRA1 protein (23-248AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine protease HTRA1 (EC 3.4.21.-) (High-temperature requirement A serine peptidase 1) (L56) (Serine protease 11), HTRA1, HTRA PRSS11
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HTRA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.