



HTRA1 Antibody, Biotin conjugated

Product Code	CSB-PA835695LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92743
Immunogen	Recombinant Human Serine protease HTRA1 protein (23-248AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine protease HTRA1 (EC 3.4.21.-) (High-temperature requirement A serine peptidase 1) (L56) (Serine protease 11), HTRA1, HTRA PRSS11
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HTRA1