



AIDA Antibody, FITC conjugated

Product Code	CSB-PA839294LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96BJ3
Immunogen	Recombinant Human Axin interactor, dorsalization-associated protein (14-306AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Axin interactor, dorsalization-associated protein (Axin interaction partner and dorsalization antagonist), AIDA, C1orf80
Immunogen Species	Homo sapiens (Human)
Research Area	Developmental Biology
Target Names	AIDA