



AIDA Antibody, FITC conjugated

Product Code	CSB-PA839294LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96BJ3
Immunogen	Recombinant Human Axin interactor, dorsalization-associated protein (14-306AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Axin interactor, dorsalization-associated protein (Axin interaction partner and dorsalization antagonist), AIDA, C1orf80
Immunogen Species	Homo sapiens (Human)
Research Area	Developmental Biology
Target Names	AIDA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.