



NECAB1 Antibody, HRP conjugated

Product Code	CSB-PA843298LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N987
Immunogen	Recombinant Human N-terminal EF-hand calcium-binding protein 1 protein (1-283AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	N-terminal EF-hand calcium-binding protein 1 (EF-hand calcium-binding protein 1) (Neuronal calcium-binding protein 1), NECAB1, EFCBP1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NECAB1