



MTMR9 Antibody, Biotin conjugated

Product Code	CSB-PA847211LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96QG7
Immunogen	Recombinant Human Myotubularin-related protein 9 protein (322-549AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myotubularin-related protein 9, MTMR9, C8orf9 MTMR8
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MTMR9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.