



NGEF Antibody, Biotin conjugated

Product Code	CSB-PA850809LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N5V2
Immunogen	Recombinant Human Ephexin-1 protein (495-710AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ephexin-1 (Eph-interacting exchange protein) (Neuronal guanine nucleotide exchange factor), NGEF
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NGEF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.