



NGEF Antibody, Biotin conjugated

Product Code	CSB-PA850809LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N5V2
Immunogen	Recombinant Human Ephexin-1 protein (495-710AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ephexin-1 (Eph-interacting exchange protein) (Neuronal guanine nucleotide exchange factor), NGEF
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NGEF