



STYX Antibody, Biotin conjugated

Product Code	CSB-PA851966LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WUJ0
Immunogen	Recombinant Human Serine/threonine/tyrosine-interacting protein (1-223AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine/tyrosine-interacting protein (Protein tyrosine phosphatase-like protein), STYX
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	STYX
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.