



RASGEF1A Antibody, FITC conjugated

Product Code	CSB-PA854076LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N9B8
Immunogen	Recombinant Human Ras-GEF domain-containing family member 1A protein (1-207AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras-GEF domain-containing family member 1A, RASGEF1A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RASGEF1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.