



Tf Antibody, FITC conjugated

Product Code	CSB-PA856363LC01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92111
Immunogen	Recombinant Mouse Serotransferrin protein (401-697AA)
Raised In	Rabbit
Species Reactivity	Mouse
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serotransferrin (Transferrin) (Beta-1 metal-binding globulin) (Siderophilin), Tf, Trf
Immunogen Species	Mus musculus (Mouse)
Research Area	Cardiovascular
Target Names	Tf
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.