



KRT12 Antibody, Biotin conjugated

Product Code	CSB-PA859510LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99456
Immunogen	Recombinant Human Keratin, type I cytoskeletal 12 protein (1-178AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Keratin, type I cytoskeletal 12 (Cytokeratin-12) (CK-12) (Keratin-12) (K12), KRT12
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KRT12
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.