



CAPN10 Antibody, Biotin conjugated

Product Code	CSB-PA864027LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HC96
Immunogen	Recombinant Human Calpain-10 protein (1-287AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calpain-10 (EC 3.4.22.-) (Calcium-activated neutral proteinase 10) (CANP 10), CAPN10, KIAA1845
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CAPN10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.