



GDF3 Antibody, FITC conjugated

Product Code	CSB-PA868295NC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NR23
Immunogen	Peptide sequence from Human Growth/differentiation factor 3 protein (202-222AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Growth/differentiation factor 3 (GDF-3), GDF3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GDF3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.