



CLSTN3 Antibody, Biotin conjugated

Product Code	CSB-PA871555LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BQT9
Immunogen	Recombinant Human Calsyntenin-3 protein (197-408AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calsyntenin-3 (Alcadein-beta) (Alc-beta), CLSTN3, CS3 KIAA0726
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	CLSTN3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.