



PDCL3 Antibody, Biotin conjugated

Product Code	CSB-PA872479LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2J4
Immunogen	Recombinant Human Phosducin-like protein 3 protein (1-239AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosducin-like protein 3 (HTPHLP) (PhPL3) (Viral IAP-associated factor 1) (VIAF-1), PDCL3, PhLP2A VIAF1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PDCL3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.