



CLDND1 Antibody, Biotin conjugated

Product Code	CSB-PA873687LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NY35
Immunogen	Recombinant Human Claudin domain-containing protein 1 protein (26-140AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Claudin domain-containing protein 1 (Membrane protein GENX-3745), CLDND1, C3orf4
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CLDND1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.