



KIF25 Antibody, Biotin conjugated

Product Code	CSB-PA883440LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UIL4
Immunogen	Recombinant Human Kinesin-like protein KIF25 protein (1-262AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kinesin-like protein KIF25 (Kinesin-like protein 3), KIF25, KNSL3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KIF25
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.