



CFHR5 Antibody, Biotin conjugated

Product Code	CSB-PA883624LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXR6
Immunogen	Recombinant Human Complement factor H-related protein 5 protein (374-569AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Complement factor H-related protein 5 (FHR-5), CFHR5, CFHL5 FHR5
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CFHR5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.