



CACNG6 Antibody, Biotin conjugated

Product Code	CSB-PA883625OD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXT2
Immunogen	Peptide sequence from Human Voltage-dependent calcium channel gamma-6 subunit protein (22-39AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Voltage-dependent calcium channel gamma-6 subunit (Neuronal voltage-gated calcium channel gamma-6 subunit), CACNG6
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CACNG6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.