



CSNK1G3 Antibody, Biotin conjugated

Product Code	CSB-PA896549LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6M4
Immunogen	Recombinant Human Casein kinase I isoform gamma-3 protein (1-204AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Casein kinase I isoform gamma-3 (CKI-gamma 3) (EC 2.7.11.1), CSNK1G3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CSNK1G3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.