



UTP11 Antibody, HRP conjugated

Product Code	CSB-PA897542LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y3A2
Immunogen	Recombinant Human Probable U3 small nucleolar RNA-associated protein 11 protein (1-253AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable U3 small nucleolar RNA-associated protein 11 (U3 snoRNA-associated protein 11) (UTP11-like protein), UTP11, UTP11L
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	UTP11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.