



NPFFR2 Antibody, Biotin conjugated

Product Code	CSB-PA015983OD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5X5
Immunogen	Peptide sequence from Human Neuropeptide FF receptor 2 protein (23-43AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neuropeptide FF receptor 2 (G-protein coupled receptor 74) (G-protein coupled receptor HLWAR77) (Neuropeptide G-protein coupled receptor), NPFFR2, GPR74 NPFF2 NPGPR
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NPFFR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.