



Trypsin Antibody, HRP conjugated

Product Code	CSB-PA018814YB01PI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00761
Immunogen	Recombinant Pig Trypsin protein (9-231AA)
Raised In	Rabbit
Species Reactivity	Pig
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Trypsin (EC 3.4.21.4)
Immunogen Species	Pig
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.