



GHRH Antibody, HRP conjugated

Product Code	CSB-PA009412EB01SH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07217
Immunogen	Recombinant Sheep Somatoliberin protein (1-44AA)
Raised In	Rabbit
Species Reactivity	Sheep
Tested Applications	ELISA
Relevance	GRF is released by the hypothalamus and acts on the adenohipophyse to stimulate the secretion of growth hormone.
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Somatoliberin (Growth hormone-releasing factor) (GRF) (Growth hormone-releasing hormone), GHRH
Species	Sheep
Research Area	Signal Transduction
Target Names	GHRH