



GHRH Antibody

Product Code	CSB-PA009412EA01SH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07217
Immunogen	Recombinant Sheep Somatoliberin protein (1-44AA)
Raised In	Rabbit
Species Reactivity	Sheep
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Somatoliberin (Growth hormone-releasing factor) (GRF) (Growth hormone-releasing hormone), GHRH
Immunogen Species	Ovis aries (Sheep)
Research Area	Signal Transduction
Target Names	GHRH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.