




# Recombinant Human Sialidase-3 (NEU3), partial

<b>Product Code</b>	CSB-EP891986HU
<b>Relevance</b>	Plays a role in modulating the ganglioside content of the lipid bilayer at the level of mbrane-bound sialyl glycoconjugates.
<b>Abbreviation</b>	Recombinant Human NEU3 protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9UQ49
<b>Alias</b>	Ganglioside sialidase;Membrane sialidaseN-acetyl-alpha-neuraminidase 3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	EEVTTCSFNSPLFRQEDDRGITYRIPALLYIPPTHTFLAFAEKRSTRRDEDALHL VLRRGLRIGQLVQWGPKPLMEATLPGHRTMNPVWEQKSGCVLFFICVR GHVTERQQIVSGRNAARLCFIYSQDAGCSWSEVRDLTEEVIKHWATFAV GPGHGIQLQSGRLVIPAYTYIYPSWFFCFQLPCKTRPHSLMIYSDDLGVTHH GRLIRPMVTVECEVAEVTGRAGHPVLYCSARTPNRCRAEALSTDHGEGFQRL ALSRQLCEPPHGCQGSVVSFRPLEIPHRCQDSSSKDAPTIQQSSPGSSLRLEE EAGTPSESWLLYSHPTSRKQRVDLGIYLNQTPLEAACWSRPWILHCGPCGYS DLAALEEEGLFGCLFECGKQECEQIAFRLFTHREILSHLQGDCTSPGRNPSQ FKSN
<b>Research Area</b>	Metabolism
<b>Source</b>	E.coli
<b>Target Names</b>	NEU3
<b>Expression Region</b>	2-428aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	52.1kDa
<b>Protein Length</b>	Partial
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.