




Recombinant Toxoplasma gondii Profilin (PRF)

Product Code	CSB-EP1510TOV
Abbreviation	Recombinant Toxoplasma gondii PRF protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q58NA1
Product Type	Recombinant Protein
Immunogen Species	Toxoplasma gondii
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	SDWDPVVKEWLVDVTGYCCAGGIANAEDGVVFAAAADDDDGWSKLYKDDHEE DTIGEDGNACGKVSINEASTIKA AVDDGSAPNGVWIGGQKYKVVVRPEKGF EYN DCTFDITMCARSKGGAHLIKTPNGSIVIALYDEEKEQDKGNSRTSALAF AEYLH QSGY
Source	E.coli
Target Names	PRF
Expression Region	2-163aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	21.4kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.