




Recombinant Human Neural cell adhesion molecule L1 (L1CAM), partial

Product Code	CSB-YP012704HU
Relevance	Cell adhesion molecule with an important role in the development of the nervous syst. Involved in neuron-neuron adhesion, neurite fasciculation, outgrowth of neurites, etc. Binds to axonin on neurons.
Abbreviation	Recombinant Human L1CAM protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P32004
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	EAIVREGGTMALSGISDFGNISATAGENYSVVSWVPKEGQC�FRFHILFKALGE EKGGASLSPQYVSYNQSSYTQWDLQPDTDYEIHLFKERMFRHQMAVKTNGT GRVRLPP
Research Area	Cell Adhesion
Source	Yeast
Target Names	L1CAM
Protein Names	Recommended name: Neural cell adhesion molecule L1-like protein Alternative name(s): Close homolog of L1 Cleaved into the following chain: 1. Processed neural cell adhesion molecule L1-like protein
Expression Region	1003-1114aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	14.5kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.