



APOC2 Antibody, Biotin conjugated

Product Code	CSB-PA001932OD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02655
Immunogen	Peptide sequence from Human Apolipoprotein C-II protein (20-38AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Apolipoprotein C-II (Apo-CII) (ApoC-II) (Apolipoprotein C2) [Cleaved into: Proapolipoprotein C-II (ProapoC-II)], APOC2, APC2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	APOC2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.