




AVIL Antibody

Product Code	CSB-PA002464LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75366
Immunogen	Recombinant Human Advillin protein (623-819AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Advillin (p92), AVIL
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	AVIL

Image



Western Blot
Positive WB detected in: HepG2 whole cell lysate
All lanes: AVIL antibody at 1:2000
Secondary
Goat polyclonal to rabbit IgG at 1/50000 dilution
Predicted band size: 93, 92 kDa
Observed band size: 93 kDa