



CALB1 Antibody, Biotin conjugated

Product Code	CSB-PA004432LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05937
Immunogen	Recombinant Human Calbindin protein (2-261AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calbindin (Calbindin D28) (D-28K) (Vitamin D-dependent calcium-binding protein, avian-type), CALB1, CAB27
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CALB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.