



GIP Antibody, HRP conjugated

Product Code	CSB-PA009434PB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09681
Immunogen	Peptide sequence from Human Gastric inhibitory polypeptide protein (113-131AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gastric inhibitory polypeptide (GIP) (Glucose-dependent insulinotropic polypeptide) (Incretin hormone), GIP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GIP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.