



GIP Antibody, HRP conjugated

Product Code	CSB-PA009434PB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09681
Immunogen	Peptide sequence from Human Gastric inhibitory polypeptide protein (113-131AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gastric inhibitory polypeptide (GIP) (Glucose-dependent insulintropic polypeptide) (Incretin hormone), GIP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GIP