



HMOX2 Antibody, Biotin conjugated

Product Code	CSB-PA010584HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30519
Immunogen	Recombinant Human Heme oxygenase 2 protein (2-293AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heme oxygenase 2 (HO-2) (EC 1.14.14.18), HMOX2, HO2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	HMOX2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.