



PRAMEF6 Antibody, Biotin conjugated

Product Code	CSB-PA018623LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5VXH4
Immunogen	Recombinant Human PRAME family member 6 protein (209-476AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PRAME family member 6, PRAMEF6
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PRAMEF6