




OPCML Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA621799LA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q14982 |
| Immunogen | Recombinant Human Opioid-binding protein/cell adhesion molecule protein (28-322AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA, WB; Recommended dilution: WB:1:1000-1:5000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Opioid-binding protein/cell adhesion molecule (OBCAM) (OPCML) (Opioid-binding cell adhesion molecule) (IgLON family member 1), OPCML, IGLON1 OBCAM |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Neuroscience |
| Target Names | OPCML |

Image



Western Blot

Positive WB detected in: U87 whole cell lysate, PC-3 whole cell lysate, HepG2 whole cell lysate, K562 whole cell lysate

All lanes: OPCML antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 39, 38, 34 kDa

Observed band size: 39 kDa