



OPCML Antibody, HRP conjugated

Product Code	CSB-PA621799LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14982
Immunogen	Recombinant Human Opioid-binding protein/cell adhesion molecule protein (28-322AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Opioid-binding protein/cell adhesion molecule (OBCAM) (OPCML) (Opioid-binding cell adhesion molecule) (IgLON family member 1), OPCML, IGLON1 OBCAM
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	OPCML
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.