



POGLUT1 Antibody, HRP conjugated

Product Code	CSB-PA818754LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NBL1
Immunogen	Recombinant Human Protein O-glucosyltransferase 1 protein (132-392AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein O-glucosyltransferase 1 (EC 2.4.1.-) (CAP10-like 46 kDa protein) (hCLP46) (KTEL motif-containing protein 1) (Myelodysplastic syndromes relative protein) (O-glucosyltransferase Rumi homolog) (hRumi) (Protein O-xylosyltransferase POGLUT1) (EC 2.4.2.26), POGLUT1, C3orf9 CLP46 KTELC1 MDSRP
Immunogen Species	Homo sapiens (Human)
Research Area	Developmental Biology
Target Names	POGLUT1