



POGLUT1 Antibody, HRP conjugated

Product Code	CSB-PA818754LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NBL1
Immunogen	Recombinant Human Protein O-glucosyltransferase 1 protein (132-392AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein O-glucosyltransferase 1 (EC 2.4.1.-) (CAP10-like 46 kDa protein) (hCLP46) (KTEL motif-containing protein 1) (Myelodysplastic syndromes relative protein) (O-glucosyltransferase Rumi homolog) (hRumi) (Protein O-xylosyltransferase POGLUT1) (EC 2.4.2.26), POGLUT1, C3orf9 CLP46 KTELC1 MDSRP
Immunogen Species	Homo sapiens (Human)
Research Area	Developmental Biology
Target Names	POGLUT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.