



POGLUT1 Antibody, Biotin conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA818754LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8NBL1 |
| Immunogen | Recombinant Human Protein O-glucosyltransferase 1 protein (132-392AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Protein O-glucosyltransferase 1 (EC 2.4.1.-) (CAP10-like 46 kDa protein) (hCLP46) (KTEL motif-containing protein 1) (Myelodysplastic syndromes relative protein) (O-glucosyltransferase Rumi homolog) (hRumi) (Protein O-xylosyltransferase POGLUT1) (EC 2.4.2.26), POGLUT1, C3orf9 CLP46 KTELC1 MDSRP |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Developmental Biology |
| Target Names | POGLUT1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |